



# ReaSLR CSM Kit - Introduction

---

## ReaMetric's Concentrated Smear Microscopy Test


For Simultaneous Sputum Liquefaction, Cell Concentration and Staining for TB screening

*“Small changes can make big difference”*

# Characteristics of an Ideal Screening Test


---

---



**Tuberculosis Control Programs need a highly Adaptable screening test with higher reliability which can be performed by minimum skilled technician without changing the current reporting method.**

# Positioning of ReaSLR CSM



**ReaSLR CSM brings Adaptability and Reliability bundled together to address the unmet needs of TB Screening test near the patient**

# Our Approach



Complete solution for TB screening using microscopy

Unmet Needs	ReaMetric's Solutions
Sputum sample processing	<ul style="list-style-type: none"><li>• ReaSLR – a ready to use sputum processing reagent</li></ul>
Sensitivity of screening test	<ul style="list-style-type: none"><li>• Concentration of MTB cells using magnetic beads</li></ul>
Reliability of reagents	<ul style="list-style-type: none"><li>• Ready to use, Unitized, and Room temperature stable reagents</li></ul>
Infrastructure	Partnering with manufacturers of Affordable and Portable instruments <ul style="list-style-type: none"><li>• LED microscope</li></ul>

These solutions are focused to enable the existing smear microscopy test to be more Reliable and Adaptable.

# Why Smear Microscopy has lower sensitivity?

- Low sample volume - Lower Sensitivity
  - Statistical error
  - Leads to higher false negative - lower interrogation volume
- Concentrated smear microscopy - Higher sample volume per smear, Enables higher sensitivity

Smear Gradation	Number of Bacilli	Estimated MTB Cells/ml	Number of HPF Interrogated	Estimated number of cells in interrogation volume	Interrogation volume at 100X or 100 fields*	Net result
Scanty	1-9	10 <sup>4</sup>	100	1	0.1µl	Statistical error rate is high; leading to more false negative
1+	10-99	10 <sup>5</sup>		10-100	0.1µl	
2+	1-10	10 <sup>6</sup>	10	1-10 (per HPF)	0.01µl	Microscopic evaluation performs well in this range
3+	>10	10 <sup>7</sup>		>100 (per HPF)	0.01µl	


AFB = Acid Fast Bacilli, HPF = High Powered Field

\* Assuming that the MTB cells are homogeneously distributed in sample. In reality, MTB cells are localized in purulent part of the sputum

Concentration of sputum sample → Higher interrogation volume → Better sensitivity

# Concentrated Smear Microscopy

**Sputum Liquefaction and Concentration is an essential for culturing the MTB**  
**Similar protocol is used for Single sputum concentrated smear microscopy**



*Workflow for concentrated smear microscopy*


**Conventional Sputum Liquefaction Method use  
Modified Petroff's Method which uses NALC-NaOH for sputum liquefaction**

**ReaMetric Introduces ReaSLR – A new Sputum Liquefaction Method  
A ready-to-use, unitized, dry reagent that is room temperature stable**

# Is ReaSLR better than NALC-NaOH?

Product Features	NALC-NaOH	ReaSLR
Liquefaction	✓	✓
Decontamination	✓	✓
Compatibility with downstream applications	✓	✓
Incubation Time (15-20 min)	✓	✓
Ready-to-use Reagent	✗	✓
Robust Process : Removal of operator / Process Variability	✗ *	✓
Process at Neutral pH	✗	✓
Processed Sample Stainability	✗	✓
Viability of MTB is maintained	✗	✓
Batch Processing	✗	✓
*prolonged incubation (>15-20 min) leads to disintegration of MTB, reducing the % viability		

# ReaSLR CSM Kit



## Kit contents: (for 100 Test)


1. **ReaSLR** - ready to use, unitized, room temperature stable sputum liquefying tablet (100 Sachets)
2. **Capture Buffer** - Magnetic nanoparticles for concentrating Mycobacteria cells (100 mL)
3. **MTB Spark** the staining solution (250 mL)
4. **Elution buffer** - For eluting out cells from TB beads (25 mL)
5. **Decolouriser solution** - To wash excess dye from the slide (250 mL)

# Training Video


---

3 min training video


# Product Performance



**AFB – ZN Staining**



**AFB – Auramine Staining**



**MTB Spark –  
Modified Auromine**

**MTB cells stained with 'MTB Spark' shows bright cells**

# Product Performance

(from ongoing vendor independent validation study)

	Direct Smear	ReaSLR CSM
Negative	76	56
Scanty	7	12
1+	3	9
2+	3	11
3+	15	16
<b>Total</b>	<b>104</b>	<b>104</b>

- Higher number of positives are detected
- 20 out of 76 Negative Direct smears are declared as ReaSLR CSM positive
- The AFB Grading improves

		Total	ReaSLR method CSM				
			Negative	Scanty	1+	2+	3+
<b>Total</b>		<b>104</b>	<b>56</b>	<b>12</b>	<b>9</b>	<b>11</b>	<b>16</b>
<b>Direct Smear</b>	Negative	76	55	9	6	4	2
	Scanty	7	1	3	2	1	0
	1+	3	0	0	1	2	0
	2+	3	0	0	0	2	1
	3+	15	0	0	0	2	13

**ReaSLR CSM method provides results in POSITIVE or NEGATIVE format (instead of AFB Grading as the sample is concentrated by 20 folds)**

**Culture result – awaiting to conclude the data with Sensitivity and Specificity**



# Thank you

[www.reametrix.com](http://www.reametrix.com)  
[enquiry@reametrix.com](mailto:enquiry@reametrix.com)